Compliance Summary for Wireless 2K Camera System

The Wireless 2K Camera System, developed for a medical device client, has undergone extensive testing and certification to ensure compliance with international regulatory standards for safety, electromagnetic compatibility, and medical device performance. Below is a detailed summary of the certifications, tests, and standards the system meets, as outlined in its Instructions for Use (IFU, Rev H).

United States Regulatory Compliance

1. FDA 510(k) Clearance

The Wireless 2K Camera System has received 510(k) clearance from the U.S. Food and Drug Administration (FDA), confirming its safety and effectiveness for use in diagnostic and operative arthroscopic procedures. The 510(k) clearance process involved demonstrating substantial equivalence to predicate devices, ensuring compliance with FDA requirements for medical devices. Updates to the IFU (Rev G, 8/23/2024) included statements to support the 510(k) submission, addressing clinical benefits and general safety and performance requirements (GSPR).

2. FCC Part 15 Compliance

The system complies with Part 15 of the Federal Communications Commission (FCC) Rules for Class B digital devices, ensuring minimal interference with radio communications.

- The system has obtained two FCC IDs for its camera head and wireless base station components, meeting FCC radiation exposure limits for an uncontrolled environment with a minimum separation distance of 20 cm between the radiator and users. It operates on specific channels (5190 MHz, 5230 MHz, 5755 MHz, 5795 MHz) and avoids restricted channels to prevent interference.
- Testing ensures the device does not cause harmful interference and can withstand interference that may cause undesired operation, as defined by 47 CFR § 2.1.

3. Electromagnetic Compatibility (EMC)

The system adheres to EMC standards for medical devices, ensuring safe operation in professional healthcare environments.

- **IEC 60601-1:** Compliance with safety requirements for fire, shock, and mechanical hazards.
- **CISPR 11:**
 - Group 1, Class A emissions, suitable for industrial and hospital environments.
 - Class A harmonic emissions per IEC 61000-3-2 and voltage fluctuations/flicker emissions per IEC 61000-3-3.
- **Electromagnetic Immunity:**
 - Electrostatic Discharge (ESD): IEC 61000-4-2 (± 8 kV contact, ± 15 kV air).

- Electrical Fast Transient/Burst: IEC 61000-4-4 (± 2 kV, 100 kHz).
 - Surge: IEC 61000-4-5 (± 1 kV line-to-line, ± 2 kV line-to-ground).
 - Voltage Dips/Interruptions: IEC 61000-4-11 (0% UT for 1 and 300 cycles, 70% UT for 30 cycles).
 - Power Frequency Magnetic Fields: IEC 61000-4-8 (30 A/m at 60 Hz).
 - Conducted RF: IEC 61000-4-6 (3 Vrms, 6 V in ISM bands, 0.15 MHz–80 MHz).
 - Radiated RF: IEC 61000-4-3 (3 V/m, 80 MHz–2.7 GHz).
 - RF Wireless Communications Immunity: Tested across multiple bands (e.g., TETRA 400, GSM 800/900, LTE, WLAN 802.11 a/n) with levels up to 28 V/m, per IEC 61000-4-3.
- The system maintains a 20 cm separation from portable RF devices to prevent performance degradation and uses shielded cables to meet Class B FCC limits.

European Union Compliance

4. CE Marking

The system bears the CE mark, indicating conformity with European Union (EU) health, safety, and environmental protection standards for medical devices. The IFU (Rev C, 2/22/2024) includes requirements for the European Economic Area (EEA), such as disposal information, clinical benefits, and GSPR compliance.

5. Waste Electrical and Electronic Equipment (WEEE) Directive

The system complies with the WEEE Directive, ensuring proper disposal of electronic components to minimize environmental impact. Users are instructed to follow local regulations for disposal, as noted in Section VI of the IFU.

International Standards for Cleaning and Sterilization

6. Cleaning and Sterilization Standards

The camera head and C-mount coupler comply with validated cleaning and sterilization processes for reusable medical devices:

- **AAMI TIR30:2011/(R)2016:** Cleaning instructions validated for reusable medical devices.
- **ANSI/AAMI ST98:2022:** Cleaning validation requirements for healthcare products.
- **ISO 22441:2022:** Low-temperature hydrogen peroxide plasma sterilization validated using the Sterrad® 100NX Standard Cycle.
- **CDC Guidelines:** Disinfection using 70-90% ethyl or isopropyl alcohol with a minimum 1-minute exposure time.
- These standards ensure the device can be safely reprocessed to prevent infection, with detailed instructions provided in Section V of the IFU.

Cybersecurity and Data Protection

7. HIPAA Compliance

The system supports compliance with the Health Insurance Portability and Accountability Act (HIPAA) through password protection for the wireless base station and secure data deletion protocols. Patient data is stored on a 1TB internal SSD with secure multi-pass deletion to prevent recovery, and USB drives must be stored securely to prevent data loss.

Additional Technical Specifications

8. Quality of Service (QoS)

The system meets stringent QoS requirements for wireless transmission:

- **Data Integrity:** Frame loss rate $\leq 0.05\%$.
- **Latency:** System latency ≤ 130 ms \leq , channel switching ≤ 220 ms.
- **Throughput:** Sustained 60 fps at 2K resolution (minimum 144 Mbps).
- **Network Accessibility:** Operates as a dedicated access point with automatic channel selection to avoid interference.
- Uses 5.1 GHz (16QAM modulation) and 5.8 GHz (5725–5850 MHz, 5150–5250 MHz) spread-spectrum communication for uncompressed video.

9. IP67 Waterproof Rating

The camera head enclosure meets the IP67 waterproof standard, ensuring it is dust-tight and protected against immersion in water up to 1 meter for 30 minutes, enhancing durability during cleaning and sterilization.

Conclusion

The Wireless 2K Camera System is designed and tested to meet rigorous U.S., EU, and international standards for medical devices. Its compliance with FDA 510(k), FCC Part 15, CE marking, EMC standards, cleaning and sterilization protocols, and HIPAA requirements ensures safety, reliability, and performance in arthroscopic procedures. The system's robust wireless communication and data management capabilities further enhance its suitability for professional healthcare environments.

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