Instant Visualization Rod: Advanced Forensic Imaging Technology

The Instant Visualization Rod revolutionizes forensic evidence collection with its non-destructive, portable, and highly sensitive imaging system. By utilizing water molecule humidification and 254nm UV light, it enhances the visibility of trace evidence like fingerprints, bloodstains, and saliva on complex surfaces, preserving DNA for further analysis. Integrated with ultra-wide-spectrum imaging (190nm–1100nm) and multi-focus fusion technology, it supports diverse evidence types across irregular surfaces using multiple light sources (UV, laser LEDs, white light). Compared to traditional methods—powder (contaminates evidence), chemical reagents (damages DNA), and laser systems (bulky and limited)—the Instant Visualization Rod offers rapid, efficient, and versatile on-site visualization, making it an essential tool for modern crime scene investigation.

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Technical Summary: Instant Visualization Rod

Aspect	Narrow: Water Molecule Humidification with Short-Wave UV	Broad: Multi-Dimensional Visualization	Comparison with Other Methods
Principle	Sprays water vapor to alter surface humidity, enhancing trace visibility (e.g., fingerprints, bloodstains) under 254nm UV light via differential adhesion and optical property changes.	Integrates water humidification, ultra-wide-spectrum imaging (190nm–1100nm), and multi-focus fusion with UV, laser LEDs, and white light for comprehensive evidence visualization.	Combines water adhesion and spectral imaging, addressing limitations of physical, chemical, and laser methods.
Advantages	- Non-Destructive: Preserves DNA using water as a clean reagent. - Fast &	- Multi-Method Integration: Visualizes diverse evidence on irregular surfaces. -	- Vs. Powder Method: No contamination, suits complex surfaces. - Vs. Chemical

Aspect	Narrow: Water Molecule Humidification with Short-Wave UV	Broad: Multi-Dimensional Visualization	Comparison with Other Methods
	Efficient: Simple spray enables rapid on-site use. - High Sensitivity: Reveals trace biological evidence on complex surfaces.	Flexible Imaging: Supports multiple evidence types with adjustable focus. - Portable: Lightweight ($\leq 1100g$), field-ready, DNA-safe.	Method: No DNA damage, simpler operation. - Vs. Laser Method: Portable, cost-effective, broader evidence scope.
Practical Examples	- Reveals aged fingerprints on a glass bottle using 254nm UV, preserving DNA. - Detects latent blood on patterned wallpaper at $\geq 1m$ distance.	- Captures palm prints on a cosmetic tube with UV and green laser, ensuring clarity on uneven surfaces. - Reveals trace saliva on fabric using UV and blue light, DNA-safe.	- Powder Method: Struggles on textured leather, risks contamination. - Chemical Method: Ninhydrin degrades DNA, slow process. - Laser Method: Fails on dark surfaces, bulky equipment. - Instant Rod: Versatile, portable, DNA-safe for all scenarios.

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